



HYGIA

Product Brochure Hygia BASE



Check In and Check Out
To help with faster Contact Tracing



E: hygia@datakrew.com
M: +65 9861 2535
W: www.hygia.care

A DATAKREW BRAND
CONFIDENTIAL DOCUMENT



HYGIA BASE – MODEL B



Completely
Touchless Device



Light and Sound
Notification



Plug and Play
Installation



0.1 Second
Response Time



Celsius and
Fahrenheit Switch



USB and Battery
Powered



HYGIA BASE – MODEL B

SPECIFICATIONS


Accuracy	± 0.2 °C (34-45 °C, place it in the operating environment for 30 minutes before use)
Automatic Alarm	Flashing + “Buzzer” sound
Measurement	Measuring distance: 5cm - 10cm
Screen	Digital display
Charging Method	USB charging or battery (18650 Li-on)
Installation Method	Nail hook, Double-sided adhesive sticking, Bracket fixing, or Screw to Tripod
Environment Temperature	10°C - 40°C (Recommended 15°C - 35°C)
Infrared Measuring Range	0°C - 50°C
Response Time	0.5 seconds
Input	5V DC
Weight	380 gram
Dimensions	170 x 115 x 140 mm
Stand By	72 hours
Counting Capacity	999

CAUTION:
Devices are not be used under direct sunlight




HYGIA BASE – MODEL B

INSTALLATION METHODS



Wall Hook



Adhesive Tape



Tripod Stand

COMPATIBILITY & CERTIFICATIONS



WARRANTY & OTHER DETAILS

- All devices comes with **1-year** warranty on product and **6 months** warranty on accessories (tripod stand, batteries etc.).
- Compatible with SafeEntry.
- Lead Time: 2 weeks.
- Transportations cost are extra.

CAUTION:

Devices are not be used under direct sunlight

CONFIDENTIAL DOCUMENT



#HYGIAAGAINSTCOVID19

ADDRESS

SINGAPORE OFFICE

Datakrew Pte. Ltd. #02-26, 1 Cleantech
Loop, Singapore 637141

INDIA OFFICE

Datakrew Pvt. Ltd. HAL 2nd Stage,
Indiranagar, Bangalore 560008